

BDL

Bridge Dental Laboratory Ltd

01

CLEAR ALIGNERS

scan@bridgedentallab.com | 01202 486113




ABOUT OUR ALIGNERS

Bridge Dental Lab provide professionally manufactured clear aligners made here in the UK. Our aim is to review each individual case with care, to align your patient's teeth safely and to offer the best possible price.

When sending impressions/scans, please attach any pre-treatment information (such as personal instructions, advice and/or images) that will aid in the preparation of the Treatment Plan. You can connect to us via Trios, iTero and carestream. All other scans must be exported as an STL/OBJ file and sent to us by email.

We keep the movement of teeth at .25mm per aligner as standard, however, in certain cases some teeth require more force. It will be highlighted in your plan when the movement of specific dentition is reaching the maximum average movement per aligner force (<.45mm)



Treatment Plan

You will receive your treatment plan within 3 working days. It will contain the information on the patients proposed Treatment Plan. This will include images, a time-lapse of each arch, and the specific instructions for any Inter-Proximal Reduction & Attachment Placement (where required), along with any necessary additional information.

Once you have reviewed the Treatment Plan, all we need is your approval either by email or phone response. You may decide that the plan provided by the laboratory could be altered. If you would like any adjustments, or the patient has cancelled the treatment, please let us know at this stage. Alterations may include but are not limited to extension/reduction in treatment duration and increase/reduction in patient tolerance. If this is the case, we ask you to state specific requests in your response to the treatment plan before approving.

On going support is available during the treatment plan process via

Email - Scan@bridgedentallab.com

Phone - 01202 486113

Whatsapp messages - 07307858772

When approval is received, we will proceed with the manufacture of the first batch of aligners (4 sets maximum, as default). Attachment templates are provided with every case that requires bonded composite attachments on certain teeth (as dictated in the Treatment Plan). You will receive the aligners within 7 working days from approval.



Revisions

As expected with clear aligner treatment, relapsing and resetting revisions are sometimes necessary to reach a desired goal. We handle such revisions on a case-by-case basis, and always consider what is best for the patient when determining the best course of action in these occasions. If revisions are required after the full treatment, we will offer 50% off per tray.



Treatment Plan

Treatment plans will expire within eight weeks if not approved. After this time, a new impression/scan is required, and a new treatment plan is to be generated should the case remain open. Cancelled or Non-approved treatment after eight weeks are subject to charges.

Please find an attached price list below.

We look forward to reviewing your next case!

50% OFF
YOUR FIRST CASE!

Get in touch for more information

Scan@bridgedentallab.com | 01202 486113