

# Connecting your AlliedStar scanner

How to send cases to Bridge Dental Lab from your AlliedStar scanner (AS 100, AS 200E, AS 260, Sensa).

AlliedStar scanners use **AS Connect** in the ScanPro software (linked to the Data Hub cloud). All AlliedStar models include free cloud storage.

## LAB CONNECTION DETAILS

Partner email: [scan@bridgedentallab.com](mailto:scan@bridgedentallab.com)

## Step-by-step

1. Create or log in to your AS Connect account at the AlliedStar portal: [allied-star.com](http://allied-star.com).
2. From the dashboard, click **My Partner**.
3. Click **Add Partner** and search for Bridge Dental Lab using [scan@bridgedentallab.com](mailto:scan@bridgedentallab.com).
4. Send the partnership request — we'll accept on our side.
5. When completing the digital Rx after a scan, choose **Bridge Dental Lab** as the recipient and submit.
6. Alternatively, export STL or PLY from the ScanPro software and email files directly to [scan@bridgedentallab.com](mailto:scan@bridgedentallab.com).

## Useful links

- AlliedStar: [allied-star.com](http://allied-star.com)

## What to include with every digital case

- **Appliance required**
- **Both arches scanned**
- **Bites** if required
- **Which practice** the case is for, if you are a multi-practice dentist

Please send any gold service scans before 11am so we can ensure your appliance will be back with you the following day.

### Need help connecting your scanner?

**Phone:** 01202 486113

**WhatsApp:** 07307 858772

**Email:** [scan@bridgedentallab.com](mailto:scan@bridgedentallab.com)

**Web:** [bridgedentallab.co.uk](http://bridgedentallab.co.uk)

**We accept scans from all major scanners with no extra fees — free 3D printed models with every appliance order.**

---

**Bridge Dental Laboratory Ltd**

Unit 10 & 11, Sea Vixen Industrial Estate  
Wilverley Road, Christchurch, Dorset BH23 3RU

01202 486113 · scan@bridgedentallab.com  
bridgedentallab.co.uk

*DAMAS-accredited · GDC-registered technicians*