

Connecting your Eletra scanner

How to send cases to Bridge Dental Lab from your Eletra Smart scanner.

Eletra scanners use **Eletra Cloud** integrated with the scanner software, plus open STL export for direct file transfer.

LAB CONNECTION DETAILS

Partner email: scan@bridgedentallab.com

Step-by-step

1. Sign in to Eletra Cloud via the scanner software, or via the URL provided in your Eletra setup pack.
2. After scanning, choose **Send to lab**.
3. If Bridge Dental Lab isn't yet listed in your partners, export the scan as an STL file instead.
4. Email the STL plus your prescription to scan@bridgedentallab.com, or upload via the Bridge Dental Lab customer portal.

Useful links

- Eletra Scanner: eletrascanner.com
- Eletra Cloud overview: eletrascanner.com/eletra-cloud

What to include with every digital case

- **Appliance required**
- **Both arches scanned**
- **Bites** if required
- **Which practice** the case is for, if you are a multi-practice dentist

Please send any gold service scans before 11am so we can ensure your appliance will be back with you the following day.

Need help connecting your scanner?

Phone: 01202 486113

WhatsApp: 07307 858772

Email: scan@bridgedentallab.com

Web: bridgedentallab.co.uk

We accept scans from all major scanners with no extra fees — free 3D printed models with every appliance order.

Bridge Dental Laboratory Ltd

Unit 10 & 11, Sea Vixen Industrial Estate
Wilverley Road, Christchurch, Dorset BH23 3RU

01202 486113 · scan@bridgedentallab.com
bridgedentallab.co.uk

DAMAS-accredited · GDC-registered technicians