

Connecting your Medit scanner

How to send cases to Bridge Dental Lab from your Medit scanner (i500, i700, i900).

Medit uses **Medit Link** as its case-sharing platform. Add Bridge Dental Lab as a partner once and you can send cases directly from the scanner software.

LAB CONNECTION DETAILS

Partner email: scan@bridgedentallab.com

Step-by-step

1. Open Medit Link on your computer or visit meditlink.com.
2. Log in with your Medit account.
3. In the left menu, click **Partners**.
4. Click **Invite Partner** (top right) and enter scan@bridgedentallab.com.
5. We'll receive and accept the partnership request.
6. When sending a case, select **Bridge Dental Lab** from the lab list and complete the prescription details.

Useful links

- Medit Link: meditlink.com
- Medit Link clinic workflow guide: it.medit.com/meditlink-workflow-for-clinics

What to include with every digital case

- **Appliance required**
- **Both arches scanned**
- **Bites** if required
- **Which practice** the case is for, if you are a multi-practice dentist

Please send any gold service scans before 11am so we can ensure your appliance will be back with you the following day.

Need help connecting your scanner?

Phone: 01202 486113

WhatsApp: 07307 858772

Email: scan@bridgedentallab.com

Web: bridgedentallab.co.uk

We accept scans from all major scanners with no extra fees – free 3D printed models with every appliance order.

Bridge Dental Laboratory Ltd
Unit 10 & 11, Sea Vixen Industrial Estate
Wilverley Road, Christchurch, Dorset BH23 3RU

01202 486113 · scan@bridgedentallab.com
bridgedentallab.co.uk
DAMAS-accredited · GDC-registered technicians